

## Reagent Preparation Sheet.

### Preparation of Phage Lysin digested sample for SeM ELISA

**Date prepared:**

**By:**

#### I. Hazard Assessment

To protect yourself from any possible hazards associated with this task wear eye protection. You should also wear latex, nitrile, or vinyl gloves and a lab coat with long sleeves. To protect your legs and feet wear closed shoes and long trousers. Do not wear sandals, shorts or a short skirt. Wash your hands before eating and when leaving the laboratory. You should review the MSDS for any chemical used in this procedure. In case of a spill with a toxic chemical remove all contaminated clothing and wash affected areas with copious quantities of water. Check location of the nearest safety shower. Eyes should be washed copiously for 15 minutes.

#### II. Reagent

Phage lysin digest of sample for SeM ELISA.

#### III. Purpose of reagent

The Phage Lysin C (PlyC) releases SeM from the wall of *Streptococcus equi* by digesting the peptidoglycan. This exposes epitopes on the SeM protein which are recognized by MAB212p (capture) and MAB211p (detector) monoclonal antibodies.

#### IV. Reconstituted vial and storage

Phage lysin was received at 5 mg/ml (i.e. 5 ug/ul) and aliquotted into 50 ul aliquots and frozen at -20 C.

#### V. To digest sample

The following method is used to digest the sample. The sample is intended to be on a cotton swab. This would be either a clinical swab, a swab placed in compost, a swab swept through the water (or over the plastic surface) of a water bucket, or a swab drawn over a blood plate to harvest streptococcal colonies

Twirl swab in 1 mL of PBS (0.15 M, pH 7.6) \_\_\_\_\_

Add 5 ul of phage lysin \_\_\_\_\_

Incubate 30 minutes at 37 C \_\_\_\_\_

Label, date and store unused digest at - 20 C \_\_\_\_\_

#### VI. Comments

